

Material List						
Type	Material Name	Manufacturer	Chemical Composition	Country of Origin	Certificate No	CE No
Metal	High Yellow Gold	Metalor	86.5% gold	Switzerland	Lot0000633699	CE1250
			11.5% platinum			
			1.4% zinc			
	Yellow Gold	Metalor	54.2% gold	Switzerland	Lot0000639949	CE1250
			31% palladium			
			4.8% argentum			
	Semi-precious	Metalor	61.5% palladium	Switzerland	Lot0000653105	CE1250
			24.5% argentum			
			2% indium			
Metal for PFM	High Yellow Gold	Argen	86% gold	U.S	Lot0000633699	CE0086
			11% platinum			
			1% Mn			
	Semi-precious	Argen	53.6% palladium	U.S	Lot0000653105	CE0086
			37.5% argentum			
			8.3% tin			
	Non-precious (Nickel-free) /Wirobond SG	Bego	<1% Ruthenium etc.	Germany	Lot13683	CE0197
			63.8% cobalt			
			24.8% chromium			
Porcelain Powder	Ceramic Powder	Creation / Vita	5.1% molybdenum	Austria / Germany	Lot8696 LOT44911	CE0483 CE0124
			silicon dioxide			
			alumina			
			sodium oxide			
			calcium oxide			

Material List						
Type	Material Name	Manufacturer	Chemical Composition	Country of Origin	Certificate No	CE No
Full Metal	Precious/Yellow Gold	Argendent	62% gold	America	Lot0000629085	CE0086
			1% palladium			
			26.5% argentum			
			8.95% copper			
	Non-precious (Nickel-free) /Wirobond SG	Bego	1.5% zinc	Germany	Lot13683	CE0197
			63.8% cobalt			
			24.8% chromium			
			5.1% molybdenum			

Material List						
Series	Material Name	Manufacturer	Chemical Composition	Country of Origin	Standard	CE No
Zirconia Crowns	Zenostar	Wieland	ZrO <sub>2</sub> +HfO <sub>2</sub> +Y <sub>2</sub> O <sub>3</sub> >99%	Germany		CE0123
			Y <sub>2</sub> O <sub>3</sub> >4,5%<6,0%			
			HfO <sub>2</sub> <5%			
			others<1%			
E-max Crowns	IPS E-max Press	Ivoclar Vivadent		Liechtenstein	LotS38422	CE0123
Procelain for Zirconia/E-max Copings	IPS E-max Ceram	Ivoclar Vivadent	lithium disilicate (LS2)	Liechtenstein	Lot20368	CE0123
			zirconium oxide (ZrO <sub>2</sub> )			

Material List						
Type	Material Name	Manufacturer	Chemical Composition	Country of Origin	Certificate No	CE No
Co-Cr Partial Framework	Non-precious (Nickel-free)/Wironit	Bego	63% cobalt	Germany	Lot13738	CE0197
			30% chromium			
			5% molybdenum			
			1% silicon			
			others<1%			
Resin Teeth	NEW ACE	YAMAHACHI		Japan	LotKD5564	CE0120
Acrylic Denture	Rapid Simplified	Vertex		Netherland	LotXT016P09	CE0120